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The purpose of the test is to evaluate whether the tested product or raw material does cause pro-sensitising effects on the involved cell model.

In the present study, we used a monocytes cell line, named THP-1, as prototypic blood-derived immunologically active cell. On these cells, we tested the expression of two costimulatory molecules, CD54 (Intercellular Adhesion Molecule) and CD86 (B7.2), using as positive control Nickel sulphate, a well-known contact sensitising agent, both in in vitro and in vivo models.

The increasing level of expression of CD54 and CD86 on monocytes is a signal of activation of the immune response derived from the exposition to a potentially sensitising contact antigen.

Functionally, the expression of co-stimulatory molecules on the dendritic cells means activation of the immunological response in terms of capability to present the antigen in the typical tissues (skin in our case), where, in vivo, the immune protective response is triggered.

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The test is carried out on a monocyte-like human line called THP-1. Cells are kept in RPMI containing 10% FCS and 2 mM glutamine.

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The sample has been eluted in phosphate buffer (1g/5ml) for 24 hours at 37°C. The eluted sample was diluted in the cell culture medium at different concentrations for the preliminary cytotoxicity screening on the cells aimed to decide the best concentration to test it without cytotoxic effects on the monocytes, in order to avoid false results. Cell medium exposed to the same experimental conditions were used as a negative control.

Following the results of the preliminary assay, the eluted sample was diluted directly in the cell medium at the two different dilutions in order to obtain the desired final concentrations in contact with the THP-1 cells in vitro. The exposure has been carried out for 24 h at 37°C with 5% CO₂.

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After the incubation with the tested substance and the controls, cells are collected, checked under the microscope for their vitality by staining with Trypan Blue dye and counting in a cell counter chamber, washed in PBS and then marked with a fluoresceinated anti- Intercellular Adhesion Molecule (CD54) or B7.2 (CD86) antibody and in both cases with a propidium iodide (PI) solution to simultaneously measure the percentage of dead cells. To this aim, a PI working solution (1 ml PI stock solution (1 mg PI/ml water) in 20 ml PI buffer (1 g glucose/l PBS without Ca/Mg)) was added to 10⁵ cells. A gate was set on the negative cells in order to exclude the dead cells. Only viable cells were included for analyses of CD54 or CD86 expression. After washing, to eliminate the excess antibody, the MFI (Mean Fluorescence Intensity) linked to the cells was evaluated by means of a flux cytometer (FACS, Fluorescence Activated Cell Sorter, Becton Dickinson, Mountain View, CA). This value is proportional to the expression of co-stimulatory molecules.

The MFI of the non-treated THP-1 cells and of cells after reaction with a monoclonal isotype-matched antibody was used as an internal control (basal fluorescence).

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The results are expressed in terms of MFI (Mean Fluorescence Intensity). MFI is the geometric mean of the fluorescence intensity of the cells decorated with the fluoresceinated antibody and it is proportional to the number of stained molecules per cell.

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Cell vitality

The percent of living cells (PI negative) must be above 80% in all samples

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For positive control

The MFI of the positive control must be > 10% of the negative control at the lowest tested concentration of the culprit substance (nichel or kathon). The MFI must display a dose-response increasing pattern with the other two tested doses.

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The costimulatory molecules expression is compared to the behaviour of Nickel, a prototypic sensitising substance, characterised by a) high increase of both the markers; b) direct correlation between concentration and intensity of the response; c) relevant effects even at very low doses.

The tested dose of 4 µg/ml of Nickel Sulphate (NiSO₄·6H₂O) corresponds to more or less 1 ppm of Nickel, dosage that is around the minimal sensitising threshold in already sensitised individuals with irritated skin. The concentration that is able to cause an allergic reaction in most of the sensitive subjects is anyway around higher value, over the 100 ppm of Nickel Sulphate·6H₂O in contact with safe and intact skin.

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The MFI values are considered significantly increased as a result of incubation with the test product, when they are higher than the value of "MFI cut-off" that is obtained by increasing the MFI of the negative control by defined percentage value.

The latter is obtained by dividing the mean + 2 standard deviations of a series of 10 determinations of MFI performed on as 10 negative control samples for the value of the MFI of the negative control itself. This determination is carried out every six months as a measure of intra-assay variability.

**Sede legale
e laboratori**Via Quarantadue Martiri 213/B 28924
Verbania (VB), ItaliaTel. +39 0323 586239
Fax +39 0323 496877
info@abich.itCF/VAT/Reg. Imp.
VCO: 01864020035
R.E.A: 189901
Cap. Soc. € 116.000 i.v.**Centro Studi Clinici
e Cosmetologici**Via Della Burrone 51, 20090
Vimodrone (MI), Italiawww.abich.it

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does not affect in this in vitro model the expression of the investigated markers in immunocompetent cells, and hence *it does not show any stimulating potential on the immune cellular response mediated by monocytes/macrophages.*

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4 **PARTE QUARTA/PART FOUR**

Bibliografia/Bibliography

- 1) EURL ECVAM RECOMMENDATION on the human Cell Line Activation Test (h-CLAT) for skin sensitisation testing. February 2015.
- 2) OECD 442E GUIDELINE FOR THE TESTING OF CHEMICALS. In vitro Skin Sensitization: human Cell Line Activation Test (h-CLAT). 29 July 2016.
- 3) E. Bocchietto, C. Paolucci, D. Breda, E. Sabbioni and SE Burastero. Human monocytoïd THP-1 cell line versus monocyte-derived human immature dendritic cells as in vitro models for predicting the sensitising potential of chemicals. International Journal of immunopathology and pharmacology Vol. 20 n°2, 261-268 (2007)
- 4) T. Ashikaga, Y Yoshida, H. Hirota et al, Development of an in vitro skin sensitization test using human cell lines: the Human Cell Line Activation Test (h-CLAT): I. Optimization of the h-CLAT protocol. Toxicology in vitro 20 (2006) 767-773.
- 5) H. Sakaguchi, T. Ashikaga, M. Miyazawa et al., Development of an in vitro skin sensitization test using human cell lines; human cell Line Activation Test (h-CLAT) II. An inter-laboratory study of the h-CLAT. Toxicology in vitro, Vol 20 Issue 5 August 2006 Pages 776-784.
- 6) Masato Mizuashi, Tomoyuki Ohtani, Satoshi Nakagawa et al, Redox imbalance induced by contact sensitizers triggers the maturation of dendritic cells. J Invest Dermatol 124:579-586, 2005.
- 7) Bocchietto E, Paolucci C, Breda D, Pecis L , Burastero S. *In vitro* prediction of contact sensitization: simultaneous evaluation of dendritic cells and of T lymphocytes in a coculture system. Atti del Congresso internazionale della Contact Dermatitis Society, Roma, 11-16 giugno 2002.
- 8) De Silva O et al. Alternative methods for skin sensitisation testing. ATLA, 1996; 24:683-705.
- 9) Thompson CB. Distinct roles for the costimulatory ligands B7-1 and B7-2 in T helper cell differentiation? Cell 1995;81:979-82.
- 10) Kuchroo VK, Das MP, Brown JA, Ranger AM, Zamvil SS, Sobel RA, Weiner HL, Nabavi N, Glimcher LH. B7-1 and B7-2 costimulatory molecules activate differentially the Th1/Th2 developmental pathways: application to autoimmune disease therapy. Cell 1995;80:707-18.
- 11) Hofer MF, Jirapongsananuruk O, Trumble AE, Leung DY. Upregulation of B7.2, but not B7.1, on B cells from patients with allergic asthma. J Allergy Clin Immunol 1998;101:96-102.
- 12) Larche M, Till SJ, Haselden BM, North J, Barkans J, Corrigan CJ, Kay AB, Robinson DS. Costimulation through CD86 is involved in airway antigen-presenting cell and T cell responses to allergen in atopic asthmatics. J Immunol 1998;161:6375-82.

- 13) Van Neerven RJ, Van de Pol MM, Van der Zee JS, Stiekema FE, De Boer M , Kapsenberg ML. Requirement of CD28-CD86 costimulation for allergen-specific T cell proliferation and cytokine expression. Clin Exp Allergy 1998;28:808-16.
- 14) Burastero SE, Magnani Z, Confetti C, Balbo P, Oddera S, Rossi GA. CD80 expression on alveolar macrophages and presence of airway eosinophils in asthma and in rhinitis before and after allergen challenge. Eur Respir J 1998;12 (suppl. 28):112s.
- 15) Burastero SE, Magnani Z, Confetti C, Abbruzzese L, Oddera S, Balbo P, Rossi GA, Crimi E. Increased expression of the CD80 accessory molecule by alveolar macrophages in asthmatic subjects and its functional involvement in allergen presentation to autologous TH2 lymphocytes. J Allergy Clin Immunol 1999;103:1136-42.
- 16) Burastero S, Magnani Z, Confetti C, Abbruzzese L, Oddera S, Balbo P, Rossi G, Crimi E. Increased expression of the CD80 accessory molecule by alveolar macrophages in asthmatic subjects and its functional involvement in allergen presentation to autologous TH2 lymphocytes. J All Clinical Immunol 1999;103:1136-1142.

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